

## Easy to use 3 button Control Unit



**- Button**  
to turn the control unit  
OFF and/or decrease  
the intensity.

**+ Button**  
to turn the control  
unit ON and/or  
increase the intensity.

**M Button**  
to select one out of  
8 available modes.

Your massage can be customized up to 19 intensities.

## Charging & General Usage

1. A red light appears when charging the control unit. The light turns off when it's fully charged.
2. Charge only when needed. Once the controller's light starts flickering, it's time to charge it. **Don't leave it plugged into an electrical outlet for long periods.**
3. If you take your feet away from the massager it will automatically turn OFF after 10-15 seconds.
4. We recommend using isopropyl alcohol with a soft cloth to clean the foot pad after use. **Please detach the control unit before cleaning the foot pad.**
5. Ensure clean, dry feet and remove excess skin or calluses before a massage for a better experience.

### **WARNING**

The neuroPro foot massager should not be used around the area of your heart if you have a pacemaker or ICD; electrical currents and electromagnetic sources may cause potential issues with the functioning of your device. Additionally, if you have deep vein thrombosis (DVT), are pregnant, or suspect you might be pregnant, you should not use the neuroPro foot massager.

# say hello to your neuroPro massager

Control Unit



Charger



USB Cable



Foot Pad



User Manual



## Quick Start Guide



Connect the circular control unit and the foot pad by fastening the two metal push connectors.



Your socks should be off before using the foot pad. **Both feet must be on the pad for it to work.**



Simply Use the "+" button to turn ON your neuroPro and then select a mode using the "M" button.



Use the "+" to increase the intensity and "-" to decrease the intensity of your massage.

Control unit automatically turns off after a 15 minute session



Press "-" on level 0 to turn OFF your massager & repeat tomorrow.



Smile you have completed your first massage session.

We recommend using your massager for 15 minutes, 1-2x a day

# Different Modes



## Mode 1: Beat

Helps neuropathy users by stimulating nerve endings in their feet.



## Mode 2: Massage

Offers a soothing & relaxing experience for neuropathy users.



## Mode 3: Activation

Designed to invigorate and wake up the feet, which is ideal for neuropathy sufferers.



## Mode 4: Training

Focuses on strengthening the foot muscles, which is beneficial for neuropathy patients.



## Mode 5: Kneading

Reduce neuropathic discomfort by targeting specific pressure points.



## Mode 6: Shaping

For neuropathy users who often contend with foot stiffness, this mode can enhance flexibility.



## Mode 7: Extrusion

Can improve circulation & relieve neuropathic symptoms like swelling & numbness



## Mode 8: Slap

Enhance blood flow, reduce numbness & provide a soothing sensation for neuropathy users


# Troubleshooting




Make sure the control unit has been charged.



Make sure the control unit is turned ON, the blue light will be illuminated when it is on & beat will illuminate on the left.



Make sure both of your feet are placed on the foot pad.



Make sure the control unit is properly attached to the foot pad. Both button type connectors should click together.

## WARNING

The neuroPro foot massager should not be used around the area of your heart if you have a pacemaker or ICD; electrical currents and electromagnetic sources may cause potential issues with the functioning of your device. Additionally, if you have deep vein thrombosis (DVT), are pregnant, or suspect you might be pregnant, you should not use the neuroPro foot massager.



Voltage: 3.7V 22mA    Temperature: 0 - 60c  
Frequency: 1-100Hz    Battery: 150mAh

### SAFE STORAGE ADVISORY

Keep away from children. Do not expose to moisture, hot places or direct sunlight. Store in a dry well ventilated area.

## Need help?


**neuroPro**



support@neuropro-us.com



help.neuropro-us.com



# neuroPro

FOOT MASSAGER

**USER GUIDE**